


Mini-invasive extracorporeal CO₂ removal system

ProLUNG[®]



ProLUNG[®] 3D
for the Aquarius[™] System



A LUNG-PROTECTIVE STRATEGY

Distributed by:

NIKKISO


ESTOR

Mini-invasive extracorporeal CO₂ removal system

ProLUNG®

A LUNG-PROTECTIVE STRATEGY

ENABLE PROTECTIVE VENTILATION



- ↓ Tidal volume/plateau pressure
- ↓ Dynamic hyperinflation
- ↓ Mechanical power
- ↓ Risk of VILI

FACILITATE EXTUBATION

- ↓ Respiratory acidosis
- ↓ Diaphragm work effort
- ↓ Inspiratory effort
- ↓ Ventilator-associated risks (VILI, VAP)





PREVENT INTUBATION

- ↓ Hypercapnia
- ↓ Respiratory acidosis
- ↓ Ventilator-associated risks (VILI, VAP)



BRIDGE TO TRANSPLANT^{10,11}

- ↓ Hypercapnia
- ↓ Respiratory acidosis
- ↓ Risk of transplantation failure



PROLUNG®

QUALITY AND INNOVATION

ProLUNG® is the reference system for mini-invasive extracorporeal CO₂ removal (ECCO₂R). ProLUNG® has all the features necessary to guarantee quality ventilatory support with a clinical rationale: high CO₂ removal capacity (VCO₂ > 100 mL/min), low invasiveness for the patient (13-14 Fr bilumen catheter).



What characteristics should the ideal ECCO₂R system have?

- ☐ High CO₂-removal performance
- ☐ Biocompatibility
- ☐ Reduced priming volume
- ☐ Prolonged kit duration
- ☐ Minimal invasiveness

Why choose ProLUNG®?



ProLUNG® 3D

- ☒ Optimal CO₂-removal capacity (VCO₂ > 100 mL/min at Q_b = 400 mL/min)
- ☒ 1.81 m² membrane in polymethylpentene (PMP) covered with phosphorylcholine
- ☒ Priming volume of 125 mL (artificial lung)
- ☒ 3 days therapy (3D)

Double lumen catheter

- ☒ Low invasiveness: double lumen catheter ≥ 13 Fr
- ☒ Femoral, jugular or subclavian vascular access



PROLUNG®


FROM THEORY TO PRACTICE

The research group of Gattinoni and Quintel at UMG carried out an animal study in 2018 to evaluate the CO₂-removal capacity of the ProLUNG® system under different conditions ¹. The study included 8 adult pigs with a body weight of 57 kg. The animals were sedated, ventilated and treated with ProLUNG® using a 13 Fr catheter. The CO₂-removal capacity of the VCO₂ ML system (membrane lung) was measured under different conditions of PaCO₂, blood flow (Qb) and medical air flow.

High CO₂-removal capacity VCO₂ > 100 mL/min

The measured VCO₂ reached a maximum value of 171 mL/min.

Typically, a removal of VCO₂ > 100 mL/min can be obtained with PaCO₂ settings between 55 and 80 mmHg, an extracorporeal blood flow (Qb) of 400 mL/min and a gas flow greater than 6 L/min.




VCO₂ (ML) as a function of blood flow (Qb) and three different settings of PaCO₂: 30, 55 e 80 mmHg.

Reduction of ventilatory burden

Minute ventilation is reduced proportionally to the quantity of CO₂ removed by ProLUNG® (VCO₂ (ML)/VCO_{2tot}).

With a PaCO₂ of 74 mmHg and pH 7.3, it was possible to remove up to 138.8 mL/min of CO₂ allowing a reduction of ventilation from 7.4 to 1.9 L/min with no complications. This corresponds to a reduction in mechanical power from 9.3 to 2.6 J/min.



Minute ventilation at 3 different settings of PaCO₂ plotted as a function of VCO₂(ML)/VCO_{2tot}.

"Minimally invasive extracorporeal CO₂ removal removes a relevant amount of CO₂ thus allowing mechanical ventilation to be significantly reduced depending on extracorporeal blood flow and inflow PCO₂. Extracorporeal CO₂ removal may provide the physiologic prerequisites for controlling ventilator-induced lung injury.

The main result of this study was that a considerable amount of CO₂ was removed by the Estor ProLUNG system using only a minimally invasive cannulation and a blood flow rate similar to that used in renal dialysis".

PROLUNG®

CLINICAL APPLICATIONS

ECCO₂R is a minimally invasive extracorporeal support for the management of ventilatory insufficiency². ECCO₂R can facilitate protective ventilation at low tidal volumes or low plateau pressures in patients in mechanical ventilation, as well as facilitating rapid extubation³. In patients undergoing non-invasive ventilation (NIV) at risk of failure, ECCO₂R can prevent the invasiveness and complications of intubation³⁻⁵.

COPD⁶⁻⁹

In COPD patients with exacerbations initially managed in NIV and at risk of failure, ProLUNG® reduces the risk of intubation, thus avoiding the associated comorbidities and the prolonged hospitalization associated with invasive mechanical ventilation. In COPD patients already in invasive mechanical ventilation, ProLUNG® contributes to protective ventilation with the aim of facilitating weaning from the ventilator.

ARDS^{6,7}

In patients with moderate ARDS where it is not possible to pursue protective ventilation because of hypercapnic respiratory acidosis, ProLUNG® allows the setting of adequate tidal volumes and plateau pressures, thus avoiding the onset of volutrauma and barotrauma.

TRANSPLANTATION^{10,11}

In all phases of lung transplantation (pre-intra-post), the use of ProLUNG® protects the lung, avoiding excessive ventilator load and allowing a better management of the transplantation procedure, thus avoiding the risk of having to resort to ECMO in an emergency.

TISSUE LESIONS¹²

In the presence of tissue lesions of the respiratory system (broncho-pleural fistulas, ruptures of the trachea or diaphragmatic lesions), the use of ProLUNG® facilitates protective ventilation.

REFRACTORY ASTHMA – EXACERBATION OF BRONCHIECTASIS^{2,13}

In patients with refractory asthma or with exacerbation of bronchiectasis, ProLUNG® facilitates protective ventilation by reducing the load induced by invasive mechanical ventilation and normalizing blood pH values.



Mini-invasive extracorporeal CO₂ removal system

ProLUNG®

A LUNG-PROTECTIVE STRATEGY

ProLUNG® PLAQUA3D

Treatment modality	Hemoperfusion
Blood flow	$Q_b \leq 450$ mL/min (Aquarius™ system)
Membrane type	Polymethylpentene covered with phosphorylcholine
Membrane surface	1,81 m ²
Priming volume	Around 250 mL (artificial lung (125 mL) + blood lines)
Sterilization	Ethylene oxide
Duration of single circuit	3 days
Rinsing and priming	2 L physiological solution with 10,000 IU of heparin
Vascular access	≥ 13 Fr double lumen central venous catheter




References


1. Duscio E et al.: Extracorporeal CO₂ Removal: The Minimally Invasive Approach, Theory, and Practice. Crit Care med. 2019.
2. Boyle AJ et al.: Extracorporeal carbon dioxide removal for lowering the risk of mechanical ventilation: research questions and clinical potential for the future. Lancet Respir Med. 2018.
3. Pisani L et al.: Management of acute hypercapnic respiratory failure. Curr Opin Crit Care. 2016.
4. Vianello A et al.: Extracorporeal CO₂ removal for refractory hypercapnia in the event of acute respiratory failure. Minerva Pneumol. 2015.
5. Vianello A et al.: Successful management of acute respiratory failure in an idiopathic pulmonary fibrosis patient using an extracorporeal CO₂ removal system. Sarcoidosis Vasc Diffuse Lung Dis. 2016.
6. Grasselli G et al.: Practical Clinical Application of an Extracorporeal Carbon Dioxide Removal System in Acute Respiratory Distress Syndrome and Acute on Chronic Respiratory Failure. ASAIO journal. 2019.
7. Hilty MP et al.: Low flow veno-venous extracorporeal CO₂ removal for acute hypercapnic respiratory failure. Minerva anesthesiologica. 2017.
8. Morelli A et al.: Extracorporeal CO₂ removal in hypercapnic patients who fail noninvasive ventilation and refuse endotracheal intubation: a case series. ESICM LIVES 2015, Berlin. 2015.
9. Pisani L et al.: Effects of Extracorporeal CO₂ Removal on Inspiratory Effort and Respiratory Pattern in Patients Who Fail Weaning from Mechanical Ventilation. Am J Respir Crit Care Med. 2015.
10. Ruberto F et al.: Low-flow veno-venous extracorporeal CO₂ removal: first clinical experience in lung transplant recipients. Int J Artif Organs. 2014.
11. Soluri-Martins A et al.: How to minimise ventilator-induced lung injury in transplanted lungs: The role of protective ventilation and other strategies. Eur J Anaesthesiol. 2015.
12. Pastore A et al.: LFVVECCO₂-R to provide "lung rest" in lesions of respiratory system: experience in one patient. SIAARTI, Rome. 2013.
13. Arcaro G et al.: The Successful Management of a Patient With Exacerbation of Non-Cystic Fibrosis Bronchiectasis and Bilateral Fibrothorax Using a Venovenous Extracorporeal CO₂ Removal System. Resp Care. 2014.


 ProLUNG 3D
EUROSETS S.r.l.
Strada Statale 12, n°143
41036 Medolla
Italy
CE 0123

 SET Waste bag &
Dialyser extension
F. M. S.p.A.
Via Farini 65 bis
13043 Cigliano (VC)
Italy
CE 0123

 ProLUNG PLAQUA3D
Assembled by Estor S.p.A.
Via Newton, 12
20016 Pero
Italy

Distributed by:
Nikkiso Belgium bv
Industriepark 6
3300 Tienen
Belgium

 Aquarius system
NIKKISO Europe GmbH
Desbrocksriede 1
30855 Langenhagen
Germany
CE 0123

 AQUALINE tubing
Haemotronic S.p.A.
Via Carreri, 16
41037 Mirandola
Italy
CE 0123

